




Assessment of EGFR agonist potency by quantitative phospho flow assay of ERK 1/2

Authors: Narasimha Simma, Shalini Kolarkar, Tej Madineni, Murali Addepalli


INTRODUCTION

EGFR activation leads to persistent induction of the classical MAPK pathway identified as Extracellular signal-regulated kinase-1 (ERK1) and ERK 2. ERK plays a fundamental role in the EGFR-driven control of epidermal proliferation and inflammation. Aberrant expression of EGFR in epidermal cells lead to hyper sensitivity & inflammation while in cancers induces proliferation. Thus, assessment of EGFR agonists and antagonists has significant role in therapeutics development.


Comparative potency determination for Test Articles, TA1, TA2, TA3 & TA4 in phosphorylation of ERK 1 & 2 in A431 cells:



Serum-starved A431 cells (epidermoid carcinoma cells having constitutive expression of EGFR) were incubated for 30 and 10 min in serum free medium at 37°C, measured the basal levels of pERK1/2 (Thr 202/ Tyr 204) by flow cytometry. Observed 85% of basal pERK1/2 in 30 min as compared to 10 min.



Representative flow cytometry histogram with an overlay of dose-dependent phosphorylation of ERK in TA 1 treated A431 cells.



Dose-response curves measuring pERK in A431 cells treated with TA1-4.

Test Article	TA1	TA2	TA3	TA4
EC ₅₀ *	12.72	3.99	0.99	1.08

*Determined by GraphPad prism ver. 9.5.0.





CONCLUSION

Application of flow cytometric analysis has identified TA3 is the most potent ($EC_{50}=0.9861$) while TA1 is the least potent ($EC_{50}= 12.72$) in phosphorylating ERK.

WHY LEXTRO

Scientists at Lextro Bio Solutions have extensive experience in working with kinase agonists and antagonists where the potencies are determined as EC_{50} and IC_{50} respectively in cell-based assays. Our cell-based kinase profile assays include both 2D as well as 3D spheroids (MCTS), Organoids, Tissue samples etc., using cell lines or patient-derived primary tumours / skin tissues.

Tools employed for kinase profiling include-

- Flow cytometry
- ELISA
- Western Blot
- Antibody arrays

Lextro Bio Solutions Pvt. Ltd.

#189, Synergy Square 1 Building, 2nd Floor, Genome Valley, Hyderabad-500 078, India
www.lextrobio.com | contact@lextrobio.com

